Abdominal Palpation & Delivery Mechanism Integrated



Product Categories: Biology Lab, Educational equipments, Paediatrics, Gynae & Obstetric Model, Reproductive Models

Product Tags: biology lab, fibre model, gynae model, paediatrics, science model **Product Page**:

https://www.labappara.com/product/abdominal-palpation-delivery-mechanism-integrated/

Product Description **Abdominal Palpation & Delivery Mechanism Integrared: Scientific Educational Model Features:**

Abdominal Palpation Demonstrate the fetus, umbilical cord and placenta of vacuum-assisted delivery, flexible fetal joint, and multiple normal and abnormal positions of fetal delivery showed;

Equipped with manual delivery mechanical parts; by manual, the whole delivery process, including engagement, descending, flexion, internal rotation, extension, restoration, external rotation and fetal delivery, can be showed;

Can practise and master the comprehensive skills of normal labor, abnormal labor (dystocia), midwifery and perineum protection;

Equipped with soft elevated cushion for Leopold Maneuver practice;

Equipped with models for cervical and vaginal examination practice before and at birth;

Component: B (Fetus for Demonstration) + C (Matrix for Manual Delivery Demonstration) + D (Cervical Changes Models) + E (Lifting "Soft Pad" for Leopold Maneuver Practice)