

## Abdominal Palpation & Delivery Mechanism Integrated



**Product Categories:** [Biology Lab](#), [Educational equipments](#), [Paediatrics](#), [Gynae & Obstetric Model](#), [Reproductive Models](#)

**Product Tags:** [biology lab](#), [fibre model](#), [gynae model](#), [paediatrics](#), [science model](#)

**Product Page:**

<https://www.labappara.com/product/abdominal-palpation-delivery-mechanism-integrated/>

### Product Description

#### **Abdominal Palpation & Delivery Mechanism Integrated: Scientific**

#### **Educational Model Features:**

Abdominal Palpation Demonstrate the fetus, umbilical cord and placenta of vacuum-assisted delivery, flexible fetal joint, and multiple normal and abnormal positions of fetal delivery showed;

Equipped with manual delivery mechanical parts; by manual , the whole delivery process, including engagement, descending, flexion, internal rotation, extension, restoration, external rotation and fetal delivery, can be showed;

Can practise and master the comprehensive skills of normal labor, abnormal labor (dystocia), midwifery and perineum protection;

Equipped with soft elevated cushion for Leopold Maneuver practice;

Equipped with models for cervical and vaginal examination practice before and at birth;

Component: B ( Fetus for Demonstration ) + C ( Matrix for Manual Delivery  
Demonstration ) +D ( Cervical Changes Models ) + E ( Lifting "Soft Pad" for Leopold  
Maneuver Practice )