

## Amber Colour Flask Volumetric (Measuring)



**Product Categories:** [Laboratory Glassware](#), [Pipettes](#), [Volumetric Flasks And Measuring Cylinder](#)

**Product Tags:** [glassware](#)

**Product Page:**

<https://www.labappara.com/product/amber-colour-flask-volumetric-measuring/>

### Product Description

Amber Color Flask Volumetric (Measuring) - Laboratory Glassware Equipment

**Amber Color Flask Volumetric (Measuring)** with interchangeable Glass Solid Stopper accuracy as per "USP (US Pharmacopoeia)" with certificate.

The flask has tight capacity tolerances and assure accurate reproducible results.

We use finest quality borosilicate glass, sourced from trusted vendors to process this flask. It is used for precise dilutions and preparation of standard solutions.

Available in different sizes and specifications, this flask is highly demanded in the market.

#### **Features:**

Optically plane bottom

PTFE stopcock

Resistant to impact

Resistant to crack

Flawless finish

Easy to clean

Resistant to chemicals

Cat. No.	Capacity ml	Tolerance ± ml.
63UA/3	5	0.02
63 UA/4	10	0.02
63 UA/5	25	0.03
63 UA/6	50	0.05
63UA/7	100	0.08
63 UA/8	200	0.10
63 UA/9	250	0.12
63 UA/10	500	0.20
63UA/11	1000	0.30