

Ammeter



Product Categories: [Lab Meters](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/ammeter/>

Product Description

Ammeter- Physics Equipments

An ammeter (from Ampere Meter) is a measuring used to measure the current in a circuit.

Electric currents are measured in amperes (A), hence the name.

Instruments used to measure smaller currents.

Dc mounted in beautiful Bakelite / acrylic stand with accuracy $\pm 0\%$ F.S.D.

With lock Type terminal packed in card box.

All common ranges available.

Available Models – MO65, MR80, MR100