

Analytical Balance



Product Categories: [Balance & Scale](#), [Chemistry Equipments](#)

Product Tags: [chemistry lab instruments](#), [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/analytical-balance/>

Product Description

Analytical Balance - Laboratory Equipment

An analytical balance (often called a "lab balance") is a class of balance designed to measure small mass in the sub-milligram range. ... Electronic analytical scales measure the force needed to counter the mass being measured rather than using actual masses.

Feature -

Aluminum alloy beam, polished wooden case, with sunmica base,

In side door case.

Nicely packed in wood case