

# Analytical funnel



**Product Categories:** [Funnels](#), [Laboratory Plastic Ware](#)  
**Product Tags:** [Laboratory Equipment](#), [Laboratory plastic ware](#)  
**Product Page:** <https://www.labappara.com/product/analytical-funnel/>

## Product Description

Analytical funnel - Scientific Laboratory Equipment

**Analytical Funnel:-** These Long Stem funnels have smooth walls, 60° cone angle and external ribs to prevent air lock. Moulded in Polypropylene, these clear funnels can also be steam autoclaved.

Material : Polypropylene

Particulars	Inside $\phi$	Stem Length	Packing
35 mm	35	37	72 Pcs
50 mm	50	50	72 Pcs
62 mm	62	67	36 Pcs
75 mm	75	77	36 Pcs
100 mm	100	95	36 Pcs