

## Animal Cages



**Product Categories:** [Animal Cage](#), [Laboratory Equipments](#), [Laboratory Plastic Ware](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#), [Manufacturer of Laboratory equipment](#), [math lab instruments](#), [physics lab instruments](#), [science lab equipment](#)

**Product Page:** <https://www.labappara.com/product/animal-cages/>

### Product Description

Animal Cages - Scientific Laboratory Equipment

Available in two different sizes, these cages provide excellent shelter for different animal species in the laboratory. The metallic top provides proper ventilation and has provisions for feeding the species. These water bottles prove to be an excellent accessory for the species kept in the cage

Material : Polypropylene / Stainless Steel

Particulars	Packing
290 X 220 X 140 mm	1 Pc
430 X 270 X 150 mm	1 Pc