

## Antibiotic Zone Reader



**Product Categories:** [Laboratory Equipments](#), [Pharmaceutical Equipments](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#)

**Product Page:** <https://www.labappara.com/product/antibiotic-zone-reader/>

### Product Description

Antibiotic Zone Reader - Laboratory Equipment

The Antibiotic zone reader provides a rapid accurate method for determining the strength of antibiotic materials by measuring the diameter of an inhibited zone in petri dish.

In use, the discs are placed on a petri-dish prepared with agar and inoculated with bacteria

The bacteria cover the entire dish except for a circular inhibition zone around, it being a function of the strength of the antibiotic

This zone reader measures the diameter of the inhibited zone to 0.02 mm. with the range of 0 to 35 mm.

Work on 230 volts A.C. 50 C/S.

It is specially used in the hospital and laboratory.

### Features:

Strong construction

Low maintenance

Rust proof body

Superior functionality