

## Aspirator Bottle



**Product Categories:** [Bottles And Assesories](#), [Laboratory Plastic Ware](#)

**Product Tags:** [chemistry lab instruments](#), [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/aspirator-bottle/>

### Product Description

Aspirator Bottle - Scientific Laboratory Equipement

Aspirator Bottle:- Aspirator Bottles are made from Polypropylene, so are much lighter than glass. Aspirator Bottles are Autoclavable, Durable and Dependable and Leak Proof. The leakproof spigot screwed to the bottle body requires only a half turn to open or close. It can be used for storing, dispensing distilled water and for other laboratory purposes. Threads should be totally disengaged before autoclaving.

**Material** : Polypropylene

#### **Features :**

Fine finish

Scratch resistance

Precise design