

## Astable Multi-vibrator (Using Transistors)



**Product Categories:** [Electronics](#), [Engineering Equipment](#), [Multivibrator & Oscillator](#), [Multivibrator & Oscillator](#)

**Product Page:**

<https://www.labappara.com/product/astable-multi-vibrator-using-transistors/>

### Product Description

Astable Multi-vibrator (Using Transistors)

The instrument comprises of the following built in parts:-

Two transistors BC 547 on front panel of the instrument.

Fixed DC regulated power supply of 12V to start 'ON' the circuit.

A combination of resistances & capacitors for biasing the circuit & taking output.