

Audio Oscillator Solid State



Product Categories: [Educational equipments](#), [Electronics](#)

Product Tags: [electronics](#)

Product Page: <https://www.labappara.com/product/audio-oscillator-solid-state/>

Product Description

Audio Oscillator Solid State - Electronics Laboratory Instrument

Features:

Accurate output

Easy operation

Superb strength

Other Details:

(Sine & Square Wave)

Audio oscillator solid state Frequency Range - 1 Hz to 100 KHz in 5 decade steps

Output Voltage 0 - 10 V p-p Continuously Variable

Selector Switch for Amplitude 0-1 / 0-10 V is provided

Selector Switch for Sine & Square Wave is provided

Output Impedance - 50 ohms