

## AUTO KARL FISCHER TITRIMETER



**Product Categories:** [Analytical Lab Instruments](#), [Spectrophotometer](#)

**Product Tags:** [Analytical Lab Equipment](#), [spectrometer](#)

**Product Page:** <https://www.labappara.com/product/auto-karl-fischer-titrimeter-3/>

### Product Description

#### **Auto Karl Fischer Titrimeter: Analytical Lab Equipment**

Buzzer Available

Timer Facility available.

Pressure filling auto zero burette with 0.05ml resolution.

Dual platinum electrode display in auto mode.

Display array of 10 LED's in line.

Type : Auto zero burette, Capacity 10 ml with glass valve

Resolution :0.05 ml

Resolution: Amber color reservoir,250 ml for longer stability. Facility to transfer reagent from burette back to reservoir.

Electrode :Twin Platinum Electrode

Vessel : Specially designed moisture tight reaction vessel. Small quantity samples can be titrated

Stirrer : Motorized Teflon Coated Magnetic Stirrer

Display : 10 segments, LED bar display to indicate level of titration

Solenoid: Special Design No. tedious solenoid valve Setting

End Point : Audio & Visual alarm

Time DELAY: 10 TO 60 SEC ADJUSTABLE

Dimension: 8 Kg approx to be arrange by user

Accuracy : +\_10 ppm(0.01 %) of 25 gm of sample