

## Autoclave (Electric)



**Product Categories:** [Autoclave & Sterilizer](#), [Autoclaves](#), [Lab Heating- Equipment](#), [Laboratory Equipments](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#), [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/autoclave-electric/>

### Product Description

#### Autoclave (Electric) - Scientific Equipment

This autoclave series is designed to cover a large field of applications for hospitals and medical centers.

Electrically operated double walled. Boiler made of stainless steel.

Stainless Steel lid having radial locking system. Working pressure steel.

Stainless Steel lid having radial locking system. Working pressure indicator, control box fitted with automatic pressure control switch, pressure gauge, steam release valve, pilot lamp and on/off switch for operation.

Made of seamless aluminum body and lid. The lid is fitted with pressure gauge and double safety valve. Fitted with 2 kW heater.

**Features:**

Prevents heat loss

Provides stable temperature

Safe & uninterrupted functioning