

Autoclave (Electrical, Vertical Single Chamber)



Product Categories: [Autoclave & Sterilizer](#), [Autoclaves](#), [Lab Heating- Equipment](#), [Laboratory Equipments](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [physics lab instruments](#)

Product Page:

<https://www.labappara.com/product/autoclave-electrical-vertical-single-chamber/>

Product Description

Autoclave (Electrical, Vertical Single Chamber) - Scientific Equipment

Available in stainless steel body.

Ring made of aluminum casting

Lid made of aluminum casting with S.S Plate.

Also available in stainless steel lid as per the customer requirement.

Available in electric/non electric (fuel based).

Also available with low water cut off device.

Also available in semi automatic (timer based).

All sizes & capacity are approximate.

Highly strong and durable vacuum and steam release valves.

Equipped with safety pressure release, weight, steam release/control valve.

Electric autoclave fitted with heating element of 2 kw single phase 220 v / 50 hz.

Good quality rubber gasket. Max. Pressure 15 psi.

Weight type safety valve and pressure regulating device for better safety and functioning.

Well designed epoxy coated outer stand.

AVAILABLE IN SIZE

Size 12" x 22"

Size 12" x 22"