

# Balance Single Pan



**Product Categories:** [Balance & Scale](#), [Chemistry Equipments](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#), [physics lab instruments](#), [science lab equipment](#)

**Product Page:** <https://www.labappara.com/product/balance-single-pan/>

## Product Description

### Balance Single Pan - Laboratory Equipment

Single pan mechanical substitution balance maintains consistent response throughout.

The useful capacity is achieved by maintaining a constant load on the balance. The subtracting mass on the same side of the beam to which the sample is added.

#### Feature -

Compact designs

Easy operation

Consistent performance

Sturdiness

Competitive rate

Low maintenance