

Basic Logic Gate



Product Categories: [Digital Electronics Lab](#), [Engineering Equipment](#)

Product Page: <https://www.labappara.com/product/basic-logic-gate/>

Product Description

Basic Logic Gate

Features :

Instrument comprises of 5V/150mA DC Regulated Power Supply for logic 1 & logic 0. Circuit diagram for OR, AND & NOT gates has been printed & components mounted on front panel to demonstrate the working of DTL gates. NAND & NOR gates can be assembled using different combinations of Gates.