

Beaker (Polypropylene)



Product Categories: [Laboratory Plastic Ware](#), [Liquid-Measurement](#)

Product Tags: [Laboratory Equipment](#), [Laboratory plastic ware](#), [Manufacturer of Laboratory equipment](#)

Product Page: <https://www.labappara.com/product/beaker-polypropylene/>

Product Description

Beaker (Polypropylene) - Scientific Laboratory Equipment

Radical Beakers, moulded in Polypropylene, have excellent clarity and very good chemical resistance. These autoclavable Beakers have prominently raised graduations on the outer surface. Pouring a Solution becomes very convenient because of its tapered spout.

Material : Polypropylene

These beakers are Moulded in Polypropylene, have excellent clarity and very good chemical resistance.

These autoclavable Beakers have prominently raised graduations on the outer surface.

Particulars	Least Count	Packing
50 ml	5	12 Pcs
100 ml	5	12 Pcs
250 ml	10	12 Pcs

500 ml	10	12 Pcs
1000 ml	25	6 Pcs
2000 ml	250	6 Pcs