

Beakers Borosilicate Glass



Product Categories: [Beakers Borosilicate Glass](#), [Laboratory Glassware](#), [Quartz Ware](#)

Product Tags: [glassware](#), [Laboratory Equipments](#)

Product Page: <https://www.labappara.com/product/beakers-borosilicate-glass/>

Product Description

Beakers Borosilicate Glass - Laboratory Glassware Equipment

A Beaker Borosilicate Glass is a simple container for stirring, mixing and heating liquids commonly used in many laboratories. Beakers, moulded in borosilicate glass. Have excellent clarity and very good chemical resistance. These Beakers have prominently raised graduations on the outer surface. Pouring a solution becomes very convenient because of its tapered spout.

Features

Fine finish

Crystal clear

Scratch resistance

Capacity - 500 ml, 800 ml, 1000 ml