

Bed Side Screen



Product Categories: [Hospital Equipment & Furniture](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [Manufacturer of Laboratory equipment](#), [math lab instruments](#), [physics lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/bed-side-screen/>

Product Description

Bed Side Screen - Hospital Equipment

Constructed with MS steel round tube having dia. 25mm.

Three/ four equal section mounted on eight 5cm castors.

Epoxy powder coated.

Also available in four fold

Overall Size - 1680H X 2440W (mm)