

## Biasing Techniques of Transistor (BJT)



**Product Categories:** [Amplifier's](#), [Amplifier's](#), [Electronics](#), [Engineering Equipment](#)

**Product Page:**

<https://www.labappara.com/product/biasing-techniques-of-transistor-bjt/>

### Product Description

Biasing Techniques of Transistor (BJT)

#### **Features :**

Instrument comprises of DC Regulated Power Supply, One NPN transistor, 2 round meters to measure the voltage and current, Different type of Resistances connected inside.