

## Binocular Research Microscope (Economy model)



**Product Categories:** [Microscope](#), [Microscopes](#)

**Product Tags:** [biology lab equipment](#)

**Product Page:**

<https://www.labappara.com/product/binocular-research-microscope-economy-mode/>

### Product Description

Binocular Research Microscope (Economy Model) - Scientific Instrument

Binocular research microscope Economy Model Robust designed, die-casted sophisticated enough to fulfill your Microscopy needs. The incomparable collection of features includes separate fine and coarse control.

A rigid and heavy stand observation Bino-Head inclined at 45° and rotatable in 360° assembled from imported prism to ensure highly accurate collimation for strain free observation.

Fixed square stage with built-in co-axial mechanical stage, graduated quadruple

revolving nose piece.

Illumination Built-in-base Halogen Lamp 6 V-20 W with electrical system and regulator for Low/Medium/High.

**Features:**

Easily adjustable lenses

Replaceable components

Easy to carry

Excellent performance

Long service life

Optical combination : 5 x, 10 x (paired) (Hug.)

Objectives 5 x, 10 x, 45 x (LS), 100 x Oil Immersion (LS) (Self Paradisaical)