

Binocular Research Microscope



Product Categories: [Microscope](#), [Microscopes](#)

Product Tags: [biology lab equipment](#)

Product Page:

<https://www.labappara.com/product/binocular-research-microscope/>

Product Description

Binocular Research Microscope - Scientific Instrument

Binocular research microscope designed, die-casted sophisticated enough to fulfill your Microscopy needs. The incomparable collection of features includes separate fine and coarse control.

A rigid and heavy stand observation Bino-Head inclined at 45° and rotatable in 360° assembled from imported prism to ensure highly accurate collimation for strain free observation.

Fixed square stage with built-in co-axial mechanical stage, graduated quadruple revolving nose piece.

Illumination Built-in-base Halogen Lamp 6 V-20 W with electrical system and regulator for Low/Medium/High.

Features:

Easily adjustable lenses

Replaceable components

Easy to carry

Excellent performance

Long service life

Optical combination : 5 x, 10 x (paired) (Hug.)

Objectives 5 x, 10 x, 45 x (LS), 100 x Oil Immersion (LS) (Self Paradiisaical)