

Bottle / Liquid Filling Machine



Product Categories: [Laboratory Equipments](#), [Pharmaceutical Equipments](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/bottle-liquid-filling-machine/>

Product Description

Bottle / Liquid Filling Machine - Laboratory Equipment

A sturdy and heavy unit made of sheet mounted on a heavy angle iron duly powder coated.

Capacity of output of water is adjustable through cam system from 50 to 300 ml.

Technical Specifications:

Filling capacity variations: 50 to 300 ml

Syringe options: Single/double

Stainless steel body

Semi automatic

Electric operated

Features:

Rust proof body

Easy to operate

Long service life