


## Burette with Pinch Cock



**Product Categories:** [Burettes Borosilicate Glass](#), [Laboratory Glassware](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#), [glassware](#)

**Product Page:** <https://www.labappara.com/product/burette-pinch-cock/>

### Product Description

Burette with Pinch Cock - Laboratory Glassware Equipment

Burette with Pinch Cock

Made of Borosilicate glass

Complete Rubber tubing

Glass jet

Pinch Clip A

#### **Features:**

Long lasting performance

Durability

Lightweight

Accurate dimensions

Sleek design

Robust construction

Very high durability

Longer functionality

## Auto-zero on Reservoir Plain with Stop-cock

Capacity	Sub. Div.
10	0.05
25	0.1
50	0.1
100	0.2