

Carboys



Product Categories: [Bottles And Asssesories](#), [Laboratory Plastic Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/carboys/>

Product Description

Carboys - Scientific Laboratory Equipment

Carboys, made of polypropylene, are much lighter than glass. These autoclavable, dependable & durable Carboys are provided with wide handles making it convenient while carrying or pouring purposes. We also offer carboy with stopcock for the esteemed clients.

Particulars Packing

10 Ltrs. 1 Pc

20 Ltrs. 1 Pc