

Centrifuge (Clinical)



Product Categories: [Centrifuges Equipment](#), [Laboratory Equipments](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/centrifuge-clinical/>

Product Description

Centrifuge (Clinical) - Scientific Laboratory Equipment

Fixed head angle rotor head with five step speed regulator, operates on 220- 240V AC, 50Hz, maximum RPM 3800. For general purpose use in pathology labs, small hospitals, schools and colleges. The centrifuge with 5 step speed regulator is provided with fixed angle rotor head.

Technical Data CH 10658/1 CH 10658/3

Max. speed 4500 RPM 4000 RPM

Max. RCF 2050 xg 1650 xg

Max. tube size 15 ml 5 ml

Max. capacity 120 ml 120 ml

Connected load 0.13 kVA 0.13 Kva

Features:

Long working life

Precisely designed

High functionality

Less maintenance

Specification:

Capacity 15 ml x 4 tube 3500 R.P.M. angle head

Capacity 15 ml x 6 tube 3500 R.P.M. angle head

Capacity 15 ml x 8 tube 3500 R.P.M. angle head