

Chemical Balance



Product Categories: [Balance & Scale](#), [Chemistry Equipments](#)

Product Tags: [chemistry lab instruments](#), [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/chemical-balance/>

Product Description

Chemical Balance - Chemistry Laboratory Equipment

Balance of great precision used in quantitative chemical analysis.

A balance consisting of a lever with two equal arms and a pan suspended from each arm.

An analytical balance is so sensitive that it can detect the mass of a single grain of a chemical

Thus, if a method of direct weighing is used, the substance ought to be added to the tared container which will hold it, directly to the pan or even to weighing paper placed on the pan.

Feature -

(Semi Analytical) with seven stone agates with sunmica base capacity 200gm.

1/5mg, Nicely packed card box.