

## Clay Pipe Triangulars



**Product Categories:** [Educational equipments](#), [Laboratory Porcelain Ware](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#)

**Product Page:** <https://www.labappara.com/product/clay-pipe-triangulars/>

### Product Description

Clay Pipe Triangulars - Scientific Laboratory Instrument

These pipe clay triangles are ideal for chemistry experiments. Can withstand high temperatures. We are a renowned name in manufacturing of Triangle Flanges. We use latest technology and comprehensive quality control systems in manufacturing our products so that our customers get exactly what they want.

SIZE:- 50 mm, 60 mm, 75 mm,

To prolong the life of porcelain products:

Follow gradual heating/cooling rates

Use an oven or hot plate and an intermediate step when quicker heating/cooling rates are required

Avoid contact of heated ware with a cold surface

Carefully inspect lab ware prior to each use and do not use any product that appears defective.