

## CD / VCD / DVD Disc Player Demonstrator



**Product Categories:** [Audio Visual Lab Demonstrators & Trainers](#), [Audio Visual Lab Demonstrators & Trainers](#), [Electronics](#), [Engineering Equipment](#)

**Product Page:**

<https://www.labappara.com/product/color-tv-demonstrator-with-14-crt/>

### Product Description

#### CD / VCD / DVD Disc Player Demonstrator

CD / VCD / DVD Player Trainer is designed to provide a live demonstrator of the working of a CD / VCD / DVD Player. The circuit / mechanical assembly is mounted on suitable board for the students to visually observe the electronic circuits and mechanical parts easily. A neat block diagram along with details of test points are provided on charts. Provision is made to measure / observe the voltages and waveforms at various strategic points through sockets.

#### **Technical Specifications:**

Power Supply : 230V AC, 50 Hz, 1phase

Functions : Play, Pause, Forward, Rewind, Stop, Open, Close, Repeat, Skip and Track select (available on Main Unit and / or Remote)

Display : Provided

Laser : Semiconductor laser

TV System : PAL

Audio Output : 2.1 Audio System

Output Socket : Audio out ; Video out

Input : USB

CD's Format : Audio, MP3 CD, VCD & DVD

Remote Control : Multi function remote

Standard Accessories : Detailed instruction manual, CD (Audio), VCD, DVD, MP3 CD, USB, A/V Cord, Mains Cord and Remote Control