

Concave Mirror



Product Categories: [Educational equipments](#), [Light & Optics](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/concave-mirror/>

Product Description

Concave Mirror - Physics Optical Equipment

A concave mirror, or converging mirror, has a reflecting surface that bulges inward (away from the incident light). Concave mirrors reflect light inward to one focal point. They are used to focus light. Mirrors are highly slivered, optically true (spray painted at the back) with a beautiful ring on its outer surface.

Diameter	Focal Length
50mm	100mm
50mm	150mm
50mm	200mm
50mm	250mm
60mm	100mm
60mm	150mm

60mm	200mm
60mm	250mm
75mm	100mm
75mm	150mm
75mm	200mm
75mm	250mm
100mm	100mm
100mm	150mm
100mm	200mm
100mm	250mm