

Convex Mirror



Product Categories: [Lenses](#), [Optical Instruments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/convex-mirror-2/>

Product Description

Convex Mirror - Scientific Laboratory Educational Equipment

A convex mirror is a spherical reflecting surface (or any reflecting surface fashioned into a portion of a sphere) in which its bulging side faces the source of light.

Automobile enthusiasts often call it a fish eye mirror while other physics texts refer to it as a diverging mirror.

Diameter	Focal length
-----------------	---------------------

50mm	10-100cm
------	----------

75mm	10-100cm
------	----------