

# Convex Mirror



**Product Categories:** [Educational equipments](#), [Light & Optics](#), [Physics Equipments](#)

**Product Tags:** [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/convex-mirror/>

## Product Description

### Convex Mirror - Physics Optical Equipment

A convex mirror is a spherical reflecting surface (or any reflecting surface fashioned into a portion of a sphere) in which its bulging side faces the source of light.

Automobile enthusiasts often call it a fish eye mirror while other physics texts refer to it as a diverging mirror.

Diameter	Focal length
50 mm	10-100 cm
75 mm	10-100 cm