

Copper Rain Guage



Product Categories: [Measurements](#), [Physics Equipments](#)

Product Tags: [Measurement](#), [Physics equipment](#)

Product Page: <https://www.labappara.com/product/copper-rain-guage/>

Product Description

Copper Rain Gauge - Physics Equipment's

Constructed from copper. It has funnel, an inner receiving vessel and an outer body all of copper, and a glass measuring cylinder. The set of rain gauge is provided with metal case having removable lid to collect rainfall and a copper calibrated rain gauge jar duly marked in centimeters.

Features:

Quality tested

Hassle free performance

High strength