

## Cork Borer (Handle)



**Product Categories:** [Chemistry Equipments](#), [Lab Supplies](#)

**Product Tags:** [chemistry lab instruments](#)

**Product Page:** <https://www.labappara.com/product/cork-borer-handle/>

### Product Description

Cork Borer (Handle) - Laboratory Instrument

A cork borer, often used in a chemistry or biology laboratory,

A metal tool for cutting a hole in a cork or rubber stopper to insert glass tubing.

Cork borers usually come in a set of nested sizes along with a solid pin for pushing the removed cork (or rubber) out of the borer.

**Available in set** - 3, 6, 12

#### **Feature -**

With Handles

Brass C.P

P German pattern.