

Crucible gooch



Product Categories: [Educational equipments](#), [Laboratory Porcelain Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/crusible-gooch/>

Product Description

Crucible gooch - Scientific Laboratory Equipment

With an aim to satisfy our clients, we are engaged in offering an extensive range of Crucible, Gooch Type. These are used in many laboratories and have superior finish. These are highly demanded in market due to their availability in various specifications.. With perforated bottom Glazed inside and outside except for outside bottom surface Without lid. We offer these products at industry leading prices to our valuable clients

Capacity :- 20ml, 25ml, 30ml, 50ml, 100ml, 200ml..

To prolong the life of porcelain products:

- Follow gradual heating/cooling rates
- Use an oven or hot plate and an intermediate step when quicker heating/cooling rates are required
- Avoid contact of heated ware with a cold surface
- Carefully inspect lab ware prior to each use and do not use any product that

appears defective