

CRYO Coders



Product Categories: [Laboratory Plastic Ware](#), [Vails](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/cryo-coders/>

Product Description

CRYO Coders - Scientific Laboratory Equipment

These Coders available in different colors can easily fit into the cap of cryo vials for color identification of different samples.

Material : Polystyrene

Features:

These Coders fit onto the top of cryogenic vial closures.

These provide quick visual identification of specific vials.

It helps in Inventory maintenance.

Available in wide range of colors.

Particulars	Packing
cryo coders	100 Pcs