

Day And Night Apparatus (With Electric)



Product Categories: [Educational equipments](#), [Geography & Working Models](#)

Product Tags: [Educational Equipments](#), [Lab Equipments](#)

Product Page: <https://www.labappara.com/product/day-night-apparatus-electric/>

Product Description

Day And Night Apparatus (With Electric) - Scientific Laboratory

Educational Equipment

An electric bulb representing sun is also fixed on the opposite side of the base.

Very interesting demo model to answer why do we have day and night.

Consisting of driven political globe with angle of inclination fitted on mica top teak wood base.

A fixed aluminum arc passing in a vertical plane around the globe is also attached to the base.