

## Deflagrating Spoon (Brass)



**Product Categories:** [Chemistry Equipments](#), [General Chemistry Labware](#)

**Product Tags:** [chemistry lab instruments](#), [science lab equipment](#)

**Product Page:** <https://www.labappara.com/product/deflagrating-spoon-brass/>

### Product Description

Deflagrating Spoon (Brass) - Laboratory Equipment

**Deflagrating-spoon (chemistry)** A long-handled spoon with a cover. Lowered into a glass vessel filled with a gas to demonstrate deflagration.

A deflagration spoon is used in chemistry experiments when the burning of [phosphorus](#), sulfur or another substance is required.

It is capable for high temperatures.

It is equipped with features to facilitate working in this environment.

Such as a long handle angled to keep heat away from the hand.

Available in:

Brass wire & Brass disc

Export quality