

DIGITAL COLORIMETER



Product Categories: [Analytical Lab Instruments](#), [Calorimeters](#)

Product Tags: [Analytical Instrument](#)

Product Page: <https://www.labappara.com/product/digital-colorimeter/>

Product Description

Digital Colorimeter - Analytical Instrument

Digital Colorimeter - Highly Stable & accurate Ideal Clinical instruments for Blood & Chemical Analysis.

Range : 400 NM to 700 NM Filters 8 High Standard Filters,

Accuracy ± 0.02 O.D.,

Reproducibility ± 0.01 O.D.,

Out Put Optical Density 0 to 1.99,

Display 2 ½ Digit LED Display, Detector Selenium Photo Cell,

Light Source 6.2 V 0.3 Amp. Tungston Filament Lamp,

Minimum Volume 1 ml solution measurement,

Power 230 V $\pm 10\%$ 50 Hz AC,

Size (LXBXH) 225 x 230 x 150 mm,

Weight 4 kg. (Approx),

Accessories Test Tubes 5 Nos, Light Source Bulb, Dust Cover, Instruction Manual.