

Digital Logic Trainer (Logic Gates, Boolean Identities & Demorgan's Theorems)



Product Categories: [Digital Electronics Lab](#), [Engineering Equipment](#)

Product Page:

<https://www.labappara.com/product/digital-logic-trainer-logic-gates-booleon-identities-demorgans-theorems/>

Product Description

Digital Logic Trainer (Logic Gates, Boolean Identities & Demorgan's Theorems)

Features :

Instrument comprises of Fixed Power supply 5V/150mA, 4 SPDT switches provided for selecting logic 1 & logic 0, 2 Red LED output indicators, circuit diagram printed for 5 OR , 5 AND & 4 NOT gates & their respective ICs placed inside the cabinet & connections brought out at sockets.