

## Digital Melting Point Apparatus



**Product Categories:** [Lab Measuring Instruments](#)

**Product Tags:** [Lab Measuring Instrument](#), [Laboratory Equipment](#)

**Product Page:**

<https://www.labappara.com/product/digital-melting-point-apparatus-2/>

### Product Description

#### **Digital Melting Point Apparatus - Laboratory Instrument**

Silicon Oil Bath, Built-in Magnetic Stirrer with electronic speed controller.

Electronic Controller for adjusting the heating rate.

Digital Display of Temperature with PT-100 Sensor.

Glare-free background light with adjustable light Intensity.

Provision for holding the melting point display Suitable to determine boiling point also and accuracy 1.0 Degree.