

## Digital Melting Point Apparatus



**Product Categories:** [Laboratory Equipments](#), [Melting Point Apparatus](#)

**Product Tags:** [chemistry lab instruments](#)

**Product Page:**

<https://www.labappara.com/product/digital-melting-point-apparatus/>

### Product Description

Digital Melting Point Apparatus - Scientific Laboratory Equipment

We have various **models** of **Digital Melting Point Apparatus** to suit the need of different customers. The Melting point **apparatus** is designed to determine the melting points with high accuracy and is supplied with a silicon oil bath with heater for heating the bath. The unit has an arrangement to hold the capillary containing the sample in uniformly heating media.