

Digital Tds Meter



Product Categories: [Analytical Lab Instruments](#)

Product Tags: [Analytical Laboratory Equipment](#)

Product Page: <https://www.labappara.com/product/digital-tds-meter-2/>

Product Description

Digital Tds Meter - Analytical Laboratory Equipment

Three digit LED display with facility

Digital TDS Meter temperature compensation and digital cell constant adjustment along with TDS cell & other accessories.

Description

TDS meter

TDS meter

Measuring range

0~1999ppm

0~19.99ppt

Accuracy

±2%FS

±2%FS

Resolution

1ppm

0.01ppt

Display

4-digital LCD with Backlit

Automatic Temperature Compensation

0°C~50°C

Environment

Temperature:0°C~50°C; Humidity:≤95%

Power Supply

4*1.5V(AG13)

Dimensions

170mm*33mm*20mm

Weight

69g