

## Dissecting Microscope



**Product Categories:** [Microscope](#), [Microscopes](#)

**Product Tags:** [biology lab equipment](#)

**Product Page:** <https://www.labappara.com/product/dissecting-microscope-5/>

### Product Description

Dissecting Microscope(DM 05) - Scientific Instrument

Dissecting Microscope with bull lens . Housed in thick plywood box.

With round base having all brass parts with 10 x and 20 x eyepieces. The stage 85 mm x 75 mm provided with spring stage clips has a glass plate with arrangement for inserting white metal plate for background. Sub stage piano-concave mirror reflects sufficient daylight & artificial light on the object. Detachable hand rests on both sides of stage provided for a study support during dissection. Complete set packed in thermocol box.

Brass parts complete with 2 eyepieces 10 x and 20 x.

#### **Optical Combination:**

Eye Piece : 10 X & 20 X Eye Piece

Object stage size : 110 mm x 115 mm

#### **Features:**

Superior functionality

Long service life

Easy maintenance

User friendly

Crystal clear lenses

Light weight