

# Distilling Apparatus Borosilicate Glass



**Product Categories:** [Distilling Apparatus Borosilicate Glass](#), [Laboratory Glassware](#)

**Product Tags:** [glassware](#), [Laboratory Equipments](#)

**Product Page:** <https://www.labappara.com/product/distilling-apparatus-borosilicate-glass/>

## Product Description

Distilling Apparatus Borosilicate Glass - Laboratory Glassware Equipment

Distilling is a process of separating the component or substances from a liquid mixture by selective evaporation .

Distillation may result in essentially complete separation (nearly pure components), or it may be a partial separation that increases the concentration of selected components of the mixture.

Consists of 500 ml.

Flask, 200 mm. Graham condenser, and stopper