


Draining Rack



Product Categories: [Clamps, Holders Racks Stands](#), [Laboratory Plastic Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/draining-rack/>

Product Description

Draining Rack - Scientific Laboratory Equipment

Draining Rack provides an excellent option for post wash draining of tubes & bottles in the laboratory. The back plate of this poly-lab Draining Rack has 20 holes into which pegs are tightened with the help of screw-nuts. Bottom of the rack is provided with a tray which collects all the waste water drained down which can further be collected in a container via drain tube. This rack not rot, stain or rust.

Material : Acrylopitrile Butadiene Styrene / Polycarbonate

Features:

Crack Resistance

Durability

Compact design

Indestructible

Particulars Packing

Draining Rack 1 Pc