

## Dropping Funnels (Cylinder)



**Product Categories:** [Funnels Borosilicate Glass](#), [Laboratory Glassware](#)

**Product Tags:** [glassware](#)

**Product Page:** <https://www.labappara.com/product/dropping-funnels-cylinder/>

### Product Description

Dropping Funnels (Cylinder) - Laboratory Glassware Equipment

Dropping Funnels Cylindrical, open top, with Screw type PTFE Needle Valve Stopcock, Stem with Cone.

A dropping funnel is a type of laboratory glassware used to transfer fluids. They are fitted with a stopcock which allows the flow to be controlled. Dropping funnels are useful for adding reagents slowly, i.e. drop-wise. This may be desirable when the quick addition of the reagent may result in side reactions, or if the reaction is too vigorous.

Capacity - 50 ml, 100 ml, 250 ml,