

E/M BY Thompson Method (Bar Magnet Method)



Product Categories: [Educational equipments](#), [Electronics](#)

Product Tags: [electronics](#)

Product Page:

<https://www.labappara.com/product/em-thompson-method-bar-magnet-method/>

Product Description

E/M BY Thompson Method (Bar Magnet Method) - Scientific Equipment

E/M BY method with build in power supply to operate cathode ray tube.

And then with one deluxe voltmeter for deflection volt .

C.R.T is mounted on teak wood stand with graph.

Another wooden stand is provided to hold bar magnets and deflection magnetometer.

E/M BY fully self controlled.

Features:-

High voltage power supply with one meter provided to measure the deflection voltage.

Complete with intensity, focus & deflection controls.

It is provided with accessories such as compass box, set of elico magnets, teak wood stand with CRT mounted, teak wood u stand with scale on both theorems for placing magnets.