

Eye Piece



Product Categories: [Educational equipments](#), [Eye Piece](#), [Light & Optics](#), [Optical Instruments](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/eye-piece/>

Product Description

Eye Piece:- Optical Instruments

Eye Piece - Physics Optical Equipment

With cross line reticulated suitable for spectrometer.

Consists of two Plano-convex lenses positioned with both convex faces toward the object.

Designed for use with low power achromatic objectives.

Our entire line of eyepieces will fit any international standard microscope with a 23 mm inner diameter eyepiece tube.