

Filter disc



Product Categories: [Educational equipments](#), [Laboratory Porcelain Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/filter-disc/>

Product Description

Filter disc - Scientific Laboratory Equipment

Features:

Latest technology

Reasonable price

Reliable

To prolong the life of porcelain products:

Follow gradual heating/cooling rates

Use an oven or hot plate and an intermediate step when quicker heating/cooling rates are required

Avoid contact of heated ware with a cold surface

Carefully inspect lab ware prior to each use and do not use any product that appears defective.